



The leader in point-of-care digital display  
mounting and mobility

THE NEXT GENERATION OF

# StyleView<sup>®</sup> Carts



Designed for Nurses.  
Engineered for IT.  
Priced for Procurement.  
Certified for Safety.

As the leader in both years of experience  
and point-of-care computer mounting innovation,  
Ergotron is uniquely qualified to assist you  
in your EMR, eMAR, PACS  
and CPOE implementation.

StyleView®

# Lean. Clean. Smart. Cart.



Other models available  
visit [ergotron.com](http://ergotron.com)



LCD Cart  
non-powered  
with autolock drawer



Notebook Cart  
66 Ah power  
with autolock drawer



Notebook Cart  
non-powered

## Features and Benefits

### Designed for Nurses

- Take it for a spin—StyleView quickly responds to caregivers' touch and effortlessly maneuvers to the point of care
- StyleView's small footprint, with ultra-smooth gliding casters, travels across crowded carpeted or hard floors with ease—perfect for space-constrained areas
- Comfortable, ergonomic data entry—work sitting or standing—helps promote caregiver wellness

- ① Effortless monitor landscape-to-portrait adjustment
- ② Twenty-five degrees of fluid independent monitor tilt
- ③ Side-to-side monitor pan improves data visibility for caregiver, patient and family
- ④ Five inches of independent monitor height adjustment, for up to 25-inch total height range—increases comfort and productivity for personalized viewing
- ⑤ Worksurface height adjustment for comfortable data entry sitting or standing with Constant Force (CF) patented motion technology
- ⑥ Lighted keyboard tray for nighttime usage
- ⑦ Ergonomic negative-tilt slide-out keyboard tray with left or right mouse tray ensures wrists remain in a neutral position during data entry—avoid repetitive-stress injuries and increase comfort and productivity
- ⑧ Electronic autolock drawer, programmable with keypad; second drawer accessory can be added to carts with a drawer
- ⑨ Four ultra-smooth gliding casters for extremely easy push/pull mobility; two locking casters ensure lock-in-place stability

### Engineered for IT

- Accommodates multiple configurations, including open notebooks, closed notebooks with monitor, CPU and monitor, thin clients, tablets and all-in-one computers
- Advanced open-architecture design supports future hardware upgrades
- Includes StyleLink™—the most advanced Power Management Software available: Allows remote configuration of entire fleet of carts and seamlessly integrates with Microsoft tools
- Integrated four-port USB hub and Ethernet ready-on function: Awaken enclosed computer from energy-saving hibernation mode with the touch of an external button
- Third-generation power system ensures robust reliability, performance and tremendous run-times—and easy to replace batteries
- Cart and power system are certified to UL 60601-1, CAN/CSA C22.2 No 601.1
- Compact true sine wave inverter power system ensures optimal performance
- Power System User Interface integrated at worksurface height (includes visual and audible battery status indicator)
- Secure, ventilated storage for notebook, thin client and CPU up to 3.75 inches (9.5 cm) thick

### Priced for Procurement

- **Lowest priced cart with most standard features**, including basket, front/rear handles, scanner holder, mouse pouch, advanced cable management, worksurface, four-port USB hub, Ethernet ready-on function and AC inlet powered coiled cord
- **Lowest total cost of ownership of any point-of-care cart.** StyleView is manufactured in an ISO certified regimen for superior quality. StyleView uses SLA batteries for long-life performance. When combined with StyleLink diagnostic tools and Ergotron's service plan, using StyleView means less maintenance and on-site technical support costs
- Extensive warranty: 5-year on cart; 2-year on power system; 1-year limited warranty on battery

### Certified for Safety

- For patient and caregiver safety, the cart and power system are certified to UL 60601-1, CAN/CSA C22.2 No 601.1 and NFPA leakage current limits
- Environmentally clean: RoHS, WEEE and EU (94/62/EC) packaging compliant
- Extremely easy to clean design promotes infection control
- Power system complies with EMI limits of FCC Part 15 Class A

© 2008 Ergotron, Inc. rev. 09/02/08 870-05-077-T Content is subject to change without notification

#### Americas Sales and Corporate Headquarters

St. Paul, MN USA  
(800) 888-8458  
+1-651-681-7600  
[www.ergotron.com](http://www.ergotron.com)  
[sales@ergotron.com](mailto:sales@ergotron.com)

#### EMEA Sales

Amersfoort, The Netherlands  
+31 33 45 45 600  
[www.ergotron.com](http://www.ergotron.com)  
[info.eu@ergotron.com](mailto:info.eu@ergotron.com)

#### APAC Sales

Singapore  
[www.ergotron.com](http://www.ergotron.com)  
[info.apac@ergotron.com](mailto:info.apac@ergotron.com)

#### Worldwide OEM Sales

[www.ergotron.com](http://www.ergotron.com)  
[info.oem@ergotron.com](mailto:info.oem@ergotron.com)

