




ERGOTRON®

IMAGINE
A BETTER READING EXPERIENCE



CARL'S TABLE®
CT15



STEVE'S STATION®
ADVANCED








STEVE'S STATION®
ESSENTIAL

Perfectly aligned monitors that adjust in unison. Sit/stand/recline worksurfaces for your precise ergonomic fit. Lighting, fan and privacy features that keep you working at your optimum all day. Only Anthro combines all of this and more into the best reading tables you can imagine.

Which table is best for you?
Use our Comparison Guide on the back.

Compare Anthro's Reading Table Families

Anthro's Reading Tables	Sit/Stand Height Range	Monitor Control	Monitor Configurations	Why You'll Love This Anthro Reading Table	Maximum Lift Capacity
 Carl's Table® CT15	30" to 49.5"	Monitor Control System™ (electric control of height and depth)	Move up to 4 monitors in unison on the Monitor Control System	A complete reading environment, Carl's Table CT15 gives you electric control to sit/stand/recline, plus monitor depth and height. Features built-in ambient lighting, acoustic barrier and exclusive White Glove Treatment.	300 lbs.
 Steve's Station® Advanced	23" to 49"	Monitor Easy Track® glides on ball bearings (manual control)	Move up to 6 monitors in unison on the Monitor Easy Track	Developed in partnership with radiologists across the U.S., Steve's Station incorporates all we've learned about keeping radiologists productive—and comfortable! Available in Single- and Dual-Surface models.	325 lbs.
 Steve's Station® Essential	22" to 54.5"	Monitor Easy Track glides on ball bearings (manual control)	Move up to 6 monitors in unison on the Monitor Easy Track	Steve's Station Essential has all the essential features you need—sit-stand adjustability and Monitor Easy Track depth control—at a very competitive price.	200 lbs.
 Monitor Easy Track®	Mounts atop your existing desk	Glides on ball bearings (manual control)	Move up to 6 monitors in unison	This stand-alone Monitor Easy Track mounts to the top of your existing desk, keeping your entire array of monitors in line and easily adjustable to your perfect focal depth.	130 lbs.
 Steve's Bench	30" to 49.5"	Optional articulating arms	Holds up to 6 monitors on articulating arms (arms not included)	At 90" wide, the Steve's Bench is built long to accommodate multiple users working side by side. Plus, it's built strong with power legs to easily lift all of your monitors and equipment.	350 lbs.

To learn more about Anthro's family of reading tables:
 call 503.853.9674 / email RadiologyPreSale@ergotron.com
 visit www.ergotron.com/radiology

Carl's Table®, Steve's Station®, Monitor Control System™ and Monitor Easy Track® are trademarks of Anthro Corporation. Anthro is a wholly owned subsidiary of Ergotron, Inc. 04/13/18

